

detect and identify



**NightOWL II**  
**LB 983 NC 100**

# NightOWL II LB 983 NC 100

Superiority in Molecular Optical Imaging

*Inspired  
by nature*



Bioluminescence imaging (BLI) and biofluorescence imaging (BFI) allow monitoring of gene expression in living organisms.



In 1989 BERTHOLD TECHNOLOGIES introduced its first low light imaging instrument – the LB 980 Luminograph. The first in-vivo gene expression experiments in plants and animals were performed on this instrument before 1993.

BLI utilizes light emitted by luciferase enzymes. Today bioluminescence markers can be tailored to any gene, enabling detailed research of gene function. BFI utilizes proteins, which fluoresce under illumination, either applied as exogenous reagents or endogenously expressed. Both BLI and BFI have contributed to the understanding of disease mechanisms and the development of new treatments.

detect and identify




# NightOWL II LB 983 NC 100

## Superiority in Molecular Optical Imaging


### Full-frame CCD cameras

Images are optically projected onto the front of a parallel array acting as the image plane. The array takes the image information and partitions the image into discrete elements. Those elements are defined by the number of pixels thus "quantizing" the image. The resulting rows of image information are then shifted in a parallel fashion to the readout register that subsequently shifts the row of information to the output as a serial stream of data. The process repeats until all rows are transferred off chip.




The image is reconstructed as dictated by the system. Since the parallel register is used for both image detection and readout, a mechanical shutter must be used to preserve image integrity.

This technology allows reliable image quantification, which is vital for comparative research.



### NightOWLcam NC 100

is an ultra sensitive back-illuminated (or backlit) CCD camera with midband coating enhancing the quantum efficiency up to 90% in the spectral range between 500 – 660 nm, which is optimal for firefly luciferase, GFP and its derivatives.




Back-illumination refers to a method of preparing the CCD sensor in a way that the photons directly strike the light-sensitive thinned back surface, in contrast to conventional CCDs where photons pass through non-light sensitive elements on the front of the CCD with a resulting loss of efficiency.

Today, these systems also have very low noise, and long exposures can therefore be used to integrate the signal over time and to obtain a usable signal.

Noise of a digital image consists of the signal noise and the camera noise which again comprises readout noise and dark noise, which is directly linked to the temperature. Efficient cooling of the array (absolute -80 to -90°C depending on the room temperature) ensures lowest light detection.


Camera	NC 100
CCD array type	back-illuminated
Grade	1
Sensitive area	13.3 x 13.3 mm
Pixel size	13 x 13 $\mu$ m
Pixel resolution	1024 x 1024
Spectral range from	300 to 1100 nm
Max. quantum efficiency	90% at 620 nm
80 % quantum efficiency	460-770 nm
Full well capacity	100.000 e-/pix
Readout noise	<3 e- rms



### Versatility and flexibility


The conditions required to image living organisms can be very different. For example, today gas anaesthesia is used for small animals but is never used in plant imaging. For plants control of light, temperature or humidity are of more importance.

In the field of infectious diseases or food processing the study of bacterial growth is the objective. In dermatology and material science the very faint luminescence from free radical oxygen species (ROS) is measured. In life science, quality control or forensic studies you need a very sensitive instrument for Western, Southern and Northern blots.



To cover all these applications BERTHOLD TECHNOLOGIES provides the very flexible low light luminescence and fluorescence imaging system NightOWL and a wide variety of accessories:


Moving of camera inside the cabinet	✓
Height correction in each position	✓
Large space inside the cabinet	✓
Easy exchange of camera	✓
Microscope and plant chamber adaption	✓
Power sockets inside the cabinet	✓
Control of interface inside the cabinet	✓
Positioning plates	✓
Macro table	✓
Flange	✓
Gas anaesthesia unit	✓
Workstation	✓
Fluorescence Reflectance Imaging	✓
Ring-light epi illumination	✓
Dual Line epi illumination	✓
Gooseneck spot illumination	✓
Transilluminators	✓
Orthogonal 3D-Imaging option	✓
Animal beds for multimodality imaging	✓



### Cabinet

The NightOWL cabinet is extremely light-tight preventing any interferences from ambient light.



NightOWL is the first imager with a motor-driven camera inside the cabinet. Optimum resolution and focus of the sample is achieved by automatic positioning of the camera according to the actual sample size.



NightOWL dark cabinet

Light-tight housing (1) containing the Peltier cooled CCD camera (2) with a motor-driven vertical adjustment of magnification (7.8), CCD chip (3), the lens (4) with a second vertical precision drive for focus adjustment (9), a fluorescence light source (5) with exchangeable filters, the sample table (6) supporting 2D and 3D objects from 35 (right drawing) to 260 mm (left drawing).

The camera can be moved from a height of 50 mm to 725 mm allowing focussing on every sample size up to 250 mm. For close-ups a macro table can be used. The camera is set up with flat field and height correction. This calibration eliminates non-uniformities caused by variations in the optical path due to height, illumination or lens effects.



# NightOWL II LB 983 NC 100


## Superiority in Molecular Optical Imaging

### Resolution

The total resolution of an image is a function of camera resolution, focal length and working distance. The software of NightOWL allows easy setting of sample size and sample height. Once the sample is defined the focus is set accordingly. The resolutions are therefore as follows, shown by some examples:

Sample size	Resolution
20 cm	200 $\mu\text{m}$
10 cm	100 $\mu\text{m}$
5 cm	50 $\mu\text{m}$
With macro table	
2 cm	20 $\mu\text{m}$
1 cm	10 $\mu\text{m}$

The closer the camera to the sample the more photons can be collected due to the spherical angle of the lens.




Sometimes small objects have to be acquired. With the macro table the magnification goes up to 5 fold. With another 5-fold digital zoom the overall magnification up to 25-fold can be possible.



The left hind leg of the anaesthetised mouse was illuminated with the gooseneck spot illumination. Only with flexible fibre optics it is possible to image the leg between macro table and camera. The distance between camera and object was only 2.5 cm (dsRed excitation at 525 nm and emission 605 nm; 3-fold digital magnification of the mouse claw).

### Additional features

The NightOWL cabinet has enough space to install special light sources or to place transilluminators, heaters, coolers etc. These devices may even be switched on and off through the software and the built-in sockets. This possibility enables the researcher to add more features into the cabinet. Plant researchers often use special lamps or flash lights in their experiments. Researchers in material science sometimes need special heating devices. The transilluminators are also connected with mains.



The flange option provides light-tight access to the inner part for tubings, cables or even fibre optics, e. g. for special illumination of plants. Such modifications of the flange can be of course customized for special purpose. BERTHOLD TECHNOLOGIES will be pleased to quote for customized flanges.



NightOWL II is equipped with a telescopic table top for easier sample handling. It is very convenient to position and check the samples outside the cabinet and then to slide them inside for acquisition.

detect and identify

### Gas anaesthesia units

During the luminescence or fluorescence image acquisition rodents have to be anaesthetized. In principle, there are two ways: intraperitoneal injection of a liquid mixture of anaesthesia (e. g. ketamine / xylazine or tribromoethanol) or anaesthesia by gaseous isoflurane.

One of the benefits of gas anaesthesia is an increased luminescent signal in rodents by a factor of two compared to tribromoethanol anaesthesia. Breathing is normal, blood pressure and ATP levels are more stable. Gas anaesthesia is less harmful, so rodents can be anaesthetized for longer periods and more often per day, which is an important advantage.

### Unit for five rodents


The TEM gas anaesthesia system has been adapted to the NightOWL. The vaporizer system works with low pressure and low flow making sure no gas is leaking from the nozzles and reducing the throughput of isoflurane. A yearly calibration of the vaporizer ensures proper functioning.



The induction box can be used for both mice and rats. In collaboration with INSERM Unité 540, Montpellier, France, a special mouse tray has been developed. Up to five mice can be anaesthetized in parallel in this tray.

The tray is temperature controlled to ensure that body temperature is kept stable during imaging. To prevent crosstalk of light emission from one rodent to the other removable barriers separate five compartments.

The anaesthesia system comes complete, but if pressured air, oxygen line or scavenging line are installed in the lab, the anaesthesia system can be modified accordingly. If any gas anaesthesia unit is already present in the lab, only the inner tray has to be ordered.



#### Order information

Complete unit for 5 rodents, 220 Volt	41930
Complete unit for 5 rodents, 110 Volt	46238
Inner tray	45941

### Single unit

The macro table for LB 983 is covered with a magnetic foil. A magnetic anaesthesia gas nozzle can be mounted in any direction for optimal mouse arrangement under the camera. This single gas nozzle is connected in the same way as the unit for 5 rodents, as is the induction box. Rats can be measured directly on the table since their body temperature is more stable.



#### Order information


Macro table with temperature control for mice	51578
Single gas nozzle for mice	53192
Single gas nozzle for rats	on request

# NightOWL II LB 983 NC 100

## Superiority in Molecular Optical Imaging

### Fluorescence Reflectance Imaging (FRI)\*


requires a range of illumination devices to excite the sample from above. The signals emitted are measured by a camera positioned above. To excite the fluorophore and measure its emission, the proper illumination and set of filters have to be chosen. BERTHOLD TECHNOLOGIES offers a complete range of filters from 340 nm up to 1100 nm.



Schematic set-up of FRI


The fluorescence under illumination is either applied as exogenous agent or will be endogenously expressed.

In biofluorescence imaging (BFI) GFP and its derivatives YFP and dsRED are used. The excitation and emission optimum of these dyes are between 470, 500 or 550 nm for excitation and 530 up to 580 nm for emission. In this spectral region melanin in skin and haemoglobin in the blood vessels absorb very strong in animals. Therefore the signal intensity will decrease rapidly the deeper the fluorescent source is in the animal.




Spectral response of Melanin and Haemoglobin.

The best spectral range for penetrating an animal is between 600 nm and 900 nm. Therefore near infrared (NIR) fluorescence is a promising technique to get better signals from deep inside the animal. Researchers are developing different dyes for this application. For example Novartis in Switzerland showed the ability to bind oxazines to beta-amyloid deposits present in Alzheimer's disease. Excitation of such dyes is done at 680 nm, emission is in the range of 720 nm.



Oxazine fluorescence in an Alzheimer disease mouse.

Another example of successful NIR fluorescence is the use of Quantum Dots® 700 or 800. The Stokes' shift of these lanthanide complexes are very high (470 nm to 700 / 800 nm), additionally lifetime of emission is long (400 ns - 400 ms).




Four different concentrations of Quantum Dots® subcutaneous injected

Another advantage of the IR-region is the low amount of autofluorescence. Other false positive signals can be detected, if too much chlorophyll (strong phosphorescent substance) is in the food. Teeth also often glow due to incorporated phosphorus, as does paper, if treated with phosphate. If plastic devices are used e. g. Petri dishes, osmotic pumps etc. a special IR cut-off filter has to be used to inhibit autofluorescence at around 800 nm.



## Optics and lamp control

The unique optical system from the LB 940 Mithras Multimode Reader has been integrated in the NightOWL II model. The light beam is kept constant for each fluorescent measurement, which is ideal with the ring-light epi illumination. If the ring-light is always set at the same height, the excitation energy on the sample will always be the same.



The lamp energy can be set by a lamp factor in the software. This allows calibration of the imaging system for each fluorophor. Comparison of the amounts of different fluorophors in one sample becomes possible.

## Fluorescence Illumination Gooseneck spot illumination

An important requirement of plant researchers is to illuminate every leaf of a plant. With the gooseneck spot illumination and its flexible fibre optics BERTHOLD TECHNOLOGIES provides an excellent tool. It is possible to bend the fibre optics in every direction.

The gooseneck spot illumination is a perfect combination together with the macro table. The space between camera and object is very small, only the gooseneck spot illumination can bring light onto the imaged area (see example page 4).

### Order information

Gooseneck spot illumination 29663

\* Some techniques for generating and/or detecting light in biological subjects are patented and may require licences from third parties. Users are advised to independently determine for themselves whether their activities infringe any valid patent.

## Ring-light epi illumination

The ring-light option for fluorescence illumination allows an even distribution of light upon the sample. The ring-light is mounted on an adjustable support stand allowing positioning at heights from 14 to 16 cm.



### Order information

Ring-light epi illumination for LB 983 51685

## Dual Line epi illumination

This epi illumination option is another alternative for fluorescence illumination to image rats or, if close-up images are required.



### Order information

Dual Line epi illumination 52295

Dual Line epi illumination  
with temperature control 53109

# NightOWL II LB 983 NC 100

## Superiority in Molecular Optical Imaging

### indiGO™ Software

The new easy-to-use indiGO™ software has been developed together with users. Well organized menus and dialogue boxes guide the user through camera set-up, image processing and image analysis.



- user defined
- quantitative analysis
- display of luminescence, fluorescence or photographic images
- contrast and image enhancement tools
- colour overlay e.g. photographic image with luminescence image, e.g. fluorescent gel with the hybridization signal, or of various fluorescent images
- line plot function
- surface plot function
- zoom function (up to 5-fold)
- definition of areas of interest and evaluation
- geometrical analysis
- arithmetic functions
- data export into spreadsheet
- raw data and processed data are filed separately (according to GLP rules)
- individual exposures or image sequences
- function to automate image processing steps
- image import and export (16-bit TIF file generated by indiGO™ can easily be processed by further software packages, e.g. for multimodality or co-registration)
- printing on any Windows printer via software
- version can be installed at instrument and office site
- remote control via internet for service and quick assistance

### Software Options


The newly developed indiGO™ software allows complete control over the hardware and offers all tools for image evaluation. There are special software needs which BERTHOLD TECHNOLOGIES offers as options, for example

- DICOM option
- 21 CRF part 11 option
- animal alignment option
- digital mouse atlas alignment option, etc.

### Multimodality

NightOWL is prepared for the sequential modus of multimodal imaging. Several animal beds for this purpose have been developed, for example combination is possible with

- MR instruments of Bruker and Philipps
- PET/SPECT instruments by ISE
- X ray-CT instruments of Scanco Medical
- Ultrasound instruments of VisualSonic



Animal bed for YAP-PET Scanner

### Multimodality software

Software packages to fuse MRI, PET, SPECT and X-ray CT data are available, since all these imaging technologies offer 3D-data. BLI and FRI data are the planar or 2D, so only the z-plane can be overlaid with the same z-plane of 3D-data. In case the orthogonal 3D-option is used, 3 z-planes can be overlaid. BERTHOLD TECHNOLOGIES will develop a solution to fuse BLI or FRI images with the other imaging technologies.

In case of the software package VINCI from Max-Planck-Institut, Cologne, fts-files can be already implemented.

## detect and identify

### Applications\*

Whole animals and plants can be imaged as well as blots, gels, microplates, cell culture dishes and arrays regardless of the luminescent or fluorescent markers used. Optical calibration ensures the comparability of all images captured with NightOWL.


Detection of weak light signals with CCD cameras can be achieved with high quantum efficiencies and extremely low noise levels to enable long exposure times. The camera and cabinet design are the key to superior imaging performance, complemented by scientific evaluation software for quantification.

Application	NC 100
Biochip	+
Bioluminescence	+++
Blot documentation	++
Chemiluminescence	+++
Colony counting	++
Fluorescence	++
Gel documentation	++
In-vivo Imaging	+++
Microplates	+
Microscopy	++
Multi-label measurements	+++

+good performance ++superior performance +++excellent performance


All applications and publications are presented on the web page: [www.berthold.com/bio](http://www.berthold.com/bio).

- In-vivo visualization of reporter gene expression in prokaryotic and eucaryotic cells, in living transgenic animals and plant.




- Visualization of bacterial growth in food
- In-vivo visualization of skin diseases in dermatology
- Research and product optimisation in varnish, paint and pigment production
- Imaging of chemiluminescence of solid polymers
- Detection of ROS (reactive oxygen species)
- Forensic Science
- Imaging of microplates: immunoassays, reporter genes detection, gene probes and phagocytosis

- In-vivo visualization of fluorophors, e.g.



Oxazines bound to beta-amyloid deposits as present in Alzheimer`s disease.

- In-vivo visualization of infectious diseases




Intraperitoneally inoculated with Salmonella enteritidis carrying a lux Operon of Xenorhabdus luminescens; exposure time: 60 sec. (Courtesy: P. Hill, Nottingham, UK).

- Study of circadian rhythms via reporter genes in living transgenic plants.



The time-course follows the rhythm of transcription from the CAB2 promoter over 48 hours.

- Gels and blots: imaging and measuring of chemiluminescent stained Southern, Northern and dot blots as well as Western blots.



\* Some techniques for generating and/or detecting light in biological subjects are patented and may require licences from third parties. Users are advised to independently determine for themselves whether their activities infringe any valid patent.

# NightOWL II LB 983 NC 100

## Technical Specification and Order information

NC 100	back illuminated, 1024 x 1024 pixel, quantum efficiency 90 % at 620 nm, sensitive from 300 to 1050 nm, dynamic range of 90 dB, cooling to absolute -80°C to -90°C depending on the room temperature.	
Resolution	Sample size	Resolution
	20 cm	200 µm
	10 cm	100 µm
	5 cm	50 µm
With macro table	2 cm	20 µm
	1 cm	10 µm
Exposure times	from 30 milliseconds to hours	
Pixel binning	variable to increase sensitivity	
Filters	4 excitation filters per slide 4 emission filters per wheel 340 nm up to 1100 nm additional filter slides/wheels available	
Light Source	75 W tungsten lamp	
Working distance	automated positioning of the camera allows working distances between 50 mm and 725 mm. For working distances below 50 mm the macro table has to be used. Connection to a microscope changes field of view also.	
Interfaces	to place transilluminators, heaters, coolers, light sources etc.	
Dimensions	122 x 60 x 40 cm (HxWxD)	
Weight	85 kg	
Regulations	CE, EN	

### Laboratory environment

Power Supply	110–240 V; 50/60 Hz; max 400 VA; minimum 3 sockets
Temperature Range	max 30°C
Humidity	10 – 80%, non condensing
PC Requirements	Pentium processor, 500 MHz (or better), CD ROM drive, 2 GB hard disk (or more), true colour 22" display, serial port, parallel port, free ethernet port (RJ-45 for service remote control), USB
Room	if gas anaesthesia is used room has to be ventilated; pressured air and scavenging line for surplus gas would be an asset.
Bench	stable to sustain 85 kg of the instrument; minimum size 120 x 50 cm (L x D)

### Order Information

### Order Number

NightOWL II LB 983 NC 100 complete incl. software	40508-30
--	----------

For accessories please see separate brochure. For further information of available filters please see filter data sheet.

BERTHOLD TECHNOLOGIES reserves the right to implement technical improvements and/or design changes without prior notice. NightOWL and indiGO are trademarks of BERTHOLD TECHNOLOGIES. Quantum Dot is a trademark of Invitrogen.



BERTHOLD TECHNOLOGIES GmbH & Co. KG

P.O. Box 100 163  
75312 Bad Wildbad  
Germany

Phone: +49 7081 177-0  
Fax: +49 7081 177-100  
E-mail: Bio@Berthold.com  
Internet: www.Berthold.com/Bio

detect and identify



## NightOWL Accessories


# NightOWL Accessories

Bioluminescence imaging (BLI) and biofluorescence imaging (BFI) allow monitoring of gene expression in living organisms\*.



The conditions required to image living organisms can be very different. For example, today gas anesthesia is used for small animals but is never used in plant imaging. For plants control of light, temperature or humidity are of more importance.

In the field of infectious diseases or food processing the study of bacterial growth is the objective. In dermatology and material science the very faint luminescence from free radical oxygen species (ROS) is measured. In life science, quality control or forensic studies you need a very sensitive instrument for Western, Southern and Northern blots.



To cover all these applications BERTHOLD TECHNOLOGIES provides the very flexible low light luminescence and fluorescence imaging system NightOWL and a wide variety of accessories:


Moving of camera inside the cabinet	✓
Height correction in each position	✓
Large space inside the cabinet	✓
Easy exchange of camera	✓
Microscope and plant chamber adaption	✓
Power sockets inside the cabinet	✓
Control of interface inside the cabinet	✓
Positioning plates	✓
Macro table	✓
Flange	✓
Gas anaesthesia unit	✓
Workstation	✓
Fluorescence Reflectance Imaging	✓
Ring-light epi illumination	✓
Dual Line epi illumination	✓
Gooseneck spot illumination	✓
Transilluminators	✓
Orthogonal 3D-Imaging option	✓
Animal beds for multimodality imaging	✓



\* Some techniques for generating and/or detecting light in biological subjects are patented and may require licences from third parties. Users are advised to independently determine for themselves whether their activities infringe any valid patent.

### Calibration plate

BERTHOLD TECHNOLOGIES offers a certified calibration plate for NightOWL, a so called secondary standard. The secondary standard requires an annual certification against a primary standard. Additionally a frame to hold neutral density filters is offered. Combinations of these filters allow linearity checks over six orders of magnitude.



#### Order information

Calibration plate	40105-10
Positioning plate for calibration plate	53613
Frame with 3 filters (OD 1, OD 2, OD 3)	53558

### Workstation

The workstation makes it very easy to move the NightOWL imaging systems to other locations. The workstation provides enough space for the anaesthesia system, computer etc. The workstation is covered with a stainless steel plate for easy cleaning.




#### Order information

Workstation	
Size: 180x60x70 cm (LxWxH)	51082

### Macro table

The closer the camera to the sample the more photons can be collected due to the spherical angle of the lens. Sometimes also small objects have to be acquired. With the macro table the magnification goes up to 5 fold. With another 5-fold digital zoom the overall magnification can be up to 25-fold.



The left hind leg of the anaesthetised mouse was illuminated with the gooseneck spot illumination. Only with flexible fibre optics it is possible to image the leg between macro table and camera. The distance between camera and object was only 2.5 cm (dsRed excitation at 525 nm and emission 605 nm; 3-fold digital magnification of the mouse claw).

The macro table for LB 983 is covered with a magnetic foil. A magnetic anaesthesia gas nozzle can be mounted in any direction for optimal animal arrangement under the camera.

#### Order information

Macro table for LB 981/LB 983	41613
Macro table with heating for LB 983	51578

### Fresnell option for micro plates

Due to parallax error imaging of microplates with normal optics is not possible. To overcome this an additional lens – a Fresnell lens – has to be added.

#### Order information

Fresnell option	53607
-----------------	-------

# NightOWL Accessories

## Bioluminescence Tomography

The only detectable signal in luminescent in-vivo imaging are the photons coming out of the object boundary, which is in practise the surface of the animal.

The simple measurement of bioluminescent signals from a mouse from different angles is not unique when constructing the set of all the solutions to this inverse problem.

Since the origin of the photons cannot be determined, a 3D reconstruction in luminescence is not possible without any assumptions.

This is the Statement of: Wang, G., Y. Li & M. Jiang (2004): Uniqueness theorems in bioluminescence tomography. Med. Phys. 31 (8): 2289 – 2299.

Furthermore, Wang, Li and Jiang mention a recently launched 3D-instrument where "a diffuse luminescent imaging tomography algorithm is used to reconstruct an internal source, coupled with a homogeneous scattering-media assumption".

They come to the conclusion, that "clearly, this approach may reveal subcutaneous depth information, but satisfactory reconstruction of a bioluminescent source distribution (both geometric and power) cannot be archived in general without compensation for the heterogeneous anatomy of the mouse".



## Orthogonal 3D-Imaging

Obviously the research community is not satisfied with 2D-imaging anymore. To have more information about depth of the signal BERTHOLD TECHNOLOGIES has developed an orthogonal 3D-Imaging option.



This option allows imaging of mice from top, left and right side. The sample is moved on the tray without changing the position. To obtain images from left or right side with the same magnification the camera position has to be lowered according to the same working distance.

With WinLight software the researcher is able to calculate roughly the volume of the illuminated spot, but it is not possible to do any reconstruction due to the inverse problem.



Correlation between luminescent signal and tumour volume during exponential proliferation. Luciferase activity and manually measured xenograft volume were plotted from five different mice that gave 19 separate measurements ( $r^2=0,8888$ ). Tumour volume ( $\text{mm}^3$ ) was estimated by the formula  $d_1 \times d_2 \times d_3 / 2$  ( $d_1$ =tumour length,  $d_2$ =width,  $d_3$ =height).

To solve the inverse problem in bioluminescence, another set of data with another technology has to be acquired. This technology may be X-ray-CT or MRI combined with optical imaging.

BERTHOLD TECHNOLOGIES has a strong focus on such multimodal imaging technologies, combined with the adequate multimodal software packages.

[Order information](#)

Orthogonal 3D-Imaging option

48005



## Positioning plates

For easy positioning of microplates or Petri dishes.




### Order information

Positioning plate for microtiter plates	41606
Positioning plate for Petri dishes (ø88mm)	41609
Positioning plate for 3D orthogonal optics	51678
Other positioning plates	on request

## Positioning plates for multimodality

Vevo 770 positioning plate allows quick transfer of the animal bed from Vevo 770 ultrasound scanner into NightOWL. Gas anaesthesia connection is possible.



### Order information

Positioning plate for VEVO 770	51674
Other positioning plates	on request

## Animal beds and holders for multimodality imaging

NMR (MRI), PET, CT, SPECT, Ultrasound and Optical Imaging are today's basic technologies applied in molecular imaging research. Each technology has its advantages and provides unique information. Researchers have to do sequential imaging or scanning in different instruments to get the required data.

To benefit from all technologies and to be able to compile and compare the information a mouse has to be kept in the same position during sequential scanning in different instruments. Therefore BERTHOLD TECHNOLOGIES developed different animal beds and holders.



Animal bed for YAP-PET Scanner




### Order information

ISE YAP-PET holder	on request
SkyScan 1178 holder	50756
Scanco vivaCT 40/70	on request
Other holders	on request

# NightOWL Accessories

## MACU® Multimodal Animal Carrier Unit

For long acquisition times – in MRI up to 4 hours – the animal has to be anaesthetized by gas and the temperature controlled. For this purpose the Medres MACU® - "Multimodal Animal Carrier Unit" - has been developed and adapted to the optical imaging system NightOWL for bioluminescence or bio-fluorescence imaging.



Multimodal Animal Carrier Unit

MACU flange mounted onto NightOWL

With MACU the animal can be imaged with different scanners in one anaesthesia session. Using a rectal temperature probe the unit is heated or cooled with water by a temperature control unit (stability < 0.2 °C). Anaesthesia is supplied by a mask in combination with the bite bar.

Surplus gas is removed with a vacuum line. ECG is non-invasively done by forepaw electrodes. For NMR investigation a 30 mm Helmholtz detector can be added to the setup and removed without interfering with the animal. The mobile operation and transportation unit supports all MACU features: body temperature, ECG, respiration, etc.

Since the NightOWL camera can be moved vertically inside the light-tight cabinet, the image size can be adjusted to the respective image taken by MRI or PET instruments.

MACU option for NightOWL consists of MACU carrier, flange and MACU plate carrier. To adapt the MACU flange, the NightOWL has to be equipped with the flange option (Order Number: 40275 for LB 981).

### Order informations

MACU carrier	47577
MACU flange	47581
MACU plate carrier	47582
MACU epi ring-light	47583

The MACU flange's diameter is 14 cm, so the carrier unit for a 12 cm MRI-bore fits easily. The MACU flange is light-tight, even though transparent tubes for water and gas are used.

On the MACU plate a guide rail is integrated for easy mounting and exact positioning of MACU. A stopper on the plate avoids crashes of the camera with MACU in case of accidental movement.

The remaining MACU setup, i. e. the tubes for hot water, the tubes connected to the gas anaesthesia system, the cables for sensors and triggers, the temperature control unit and the plug & work master box with low-noise power supply, is identical, regardless, if it is used in different MRI or PET instruments or the NightOWL.

### Options in combination with other imaging technologies

Optical sensor	47578
Platin-electrodes	47579
Temperature control unit	47580
Plug & Work masterbox	47567
MACU plug-in	47568
Respiration amplifier plug-in	47569
Electrocardiogram amplifier plug-in	47570
Bloodpressure amplifier plug-in	47571
Electrocardiogram trigger plug-in	47572
Respiration trigger plug-in	47573
Difference amplifier plug-in	47574
Stimulation plug-in	47575



MACU-adapter for human MRI scanner.

## Multimodality software

VINCI ("Volume Imaging in Neurological Research, Co-Registration and ROIs included") software package, developed by Max-Planck institute for neurological research, Cologne (Germany), is highly modular, expandable, compact, entirely true colour based and allows online fusion and contour rendering of several images, more than 50 studies can be displayed simultaneously in orthogonal views on current PCs. Since VINCI version 2.05, fts-files generated by NightOWL can be visualized and analyzed further.

The Swiss company PMOD ([www.pmod.com](http://www.pmod.com)) offers another software package for image registration and fusion.

## detect and identify

### Sterile boxes

Sterile boxes, offered by the company Summit can be used in NightOWL. Different options are available to fit the NightOWL flange. For small sterile box BERTHOLD TECHNOLOGIES offers a positioning plate.



#### Order information

Big sterile box	on request
Small sterile box	on request
Positioning plate for small sterile box	52638
Gas connectors for flange	47969

### Filters

BERTHOLD TECHNOLOGIES offers a complete range of filters between 340nm up to 1100nm. Filters have to be used in fluorescence and as well BRET or FRET applications. If a transilluminator is used only an emission filter is needed.


BERTHOLD TECHNOLOGIES offers filter pairs for NightOWL for following standard dyes amongst others:

	Excitation	Order info	Emission	Order info
GFP	475/20	53183	520/10	39805
Cy5	630/20	50097	680/30	49180
Cy5.5	630/20	50097	700/20	50479
Cy7	700/20	50475	780/20	50476
ICG	740/30	50480	820/30	50481
dsRed	530/20	38536	600/20	50477
Qdot700®	630/20	50097	700/20	50479
Qdot800®	630/20	50097	820/30	50481
mCherry	550/10	39796	620/10	40540

### Transilluminators

For the NightOWL system BERTHOLD TECHNOLOGIES offers special transilluminators with illuminated area on the left side ensuring a complete field of view as the camera is also positioned on the left side.

The area of all transilluminators for NightOWL is 20x20 cm. The illumination by eight lamps each with 8 Watts gives the most uniform distribution of light. Housings of stainless steel and glass cover sealed with silicon ensure stability and long lifetime.



Since NightOWL cameras are very sensitive BERTHOLD TECHNOLOGIES offers only transilluminators with adjustable intensity in 10 % steps.


#### Available transilluminators

UV transilluminator, 254 nm, intensity adjustable,	230 V	42600
	150 V	45133
UV transilluminator, 312 nm, intensity adjustable,	230 V	42601
	115 V	45134
UV transilluminator, 365 nm, intensity adjustable,	230 V	42602
	150 V	45135
Blue light transilluminator, 470 nm, intensity adjustable,	230 V	42604
	150 V	45136
Yellow light transilluminator, 595 nm, intensity adjustable,	230 V	50973
	150 V	50982
Red light transilluminator, 625 nm, intensity adjustable,	230 V	49742
	150 V	50981
White light transilluminator, intensity adjustable,	230 V	42606
	150 V	45394
Converter plates		on request


# NightOWL Accessories

## BERTHOLD TECHNOLOGIES, one of the pioneers in Molecular Imaging

BERTHOLD TECHNOLOGIES GmbH & Co. KG is located in Bad Wildbad, Germany. The company was founded in 1949 by Prof. Dr. Rudolf Berthold and was named "Laboratorium Prof. Dr. Rudolf Berthold". At the end of the 70s BERTHOLD developed an animal based imager based on a TLC scanner.



Scan of a lung of a rat using 150 µCi <sup>99m</sup>Tc with Ironhydroxid.



In the early 80s Siemens Germany, developed a high sensitive intensified camera. Prof. Szalay pioneered the non-invasive measurement of gene expression in living organisms using this camera.

In collaboration with Siemens BERTHOLD introduced the first low light imaging instrument for this application in 1989, the LB 980 Luminograph. The first in-vivo gene expression experiments in plants and animals performed on this instrument date back to the year 1993.

With the improvement of slow scanned CCD cameras and the inherent problems of linearity and dynamic range of intensified cameras BERTHOLD introduced the NightOWL in 1996. This instrument used the state of the art slow scan cooled CCD technology to provide an extremely sensitive molecular imaging system.

Today a new generation of slow scan CCD cameras have been developed further enhancing the sensitivity. Furthermore, based on the microplate reader technology of BERTHOLD, in LB 983 NightOWL II fluorescence imaging has been optimised using a sophisticated lamp and beam control with automated filter changing.

BERTHOLD TECHNOLOGIES reserves the right to implement technical improvements and/or design changes without prior notice. NightOWL is a trademark of BERTHOLD TECHNOLOGIES, Image-Pro of MediaCybernetics. MACU® is a registered trademark of medres. Quantum Dots® is a registered trademark of Invitrogen. Third party products may not be available in all countries.



BERTHOLD TECHNOLOGIES GmbH & Co. KG

P.O. Box 100 163  
75312 Bad Wildbad  
Germany

Phone: +49 7081 177-0  
Fax: +49 7081 177-100  
E-mail: [Bio@Berthold.com](mailto:Bio@Berthold.com)  
Internet: [www.Berthold.com/Bio](http://www.Berthold.com/Bio)