## COAGULATION - COAGULACIÓN

### *Clot - SP*

Versatile **coagulometer controlled by microprocessor** with turbodensidometric detection of the beginning of the Clot formation (fibrine polymerization).

The clot detection is carried out in positive O.D. variations or even in negative ones in high level fibrinogenemias.

It is very easy to use and can store up to 10 methods in memory.

Performance of all the commonly used tests like the Quick test (PT), aPTT, TT, Fibrinogen, ATIII, factor assays, etc, and so on, with a high level of precision and repeatability.

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## *Clot - SP*

### ONE CHANNEL COAGULOMETER

The detection lamp intensity is automatically adjusted after addition of the initial reagent on the measuring cell. Consequently the baseline is adjusted as well.

Automatic Clot detection is carried out in a homogeneous mixed medium by means of a system that does not interfere with the normal Clot formation and avoids possible sedimentation phenomena of activating reagents during the test, with the turbodetection concentrated in the central axis of the cuvette.

All results are displayed on the display and automatically printed out on continuous paper, together with the ID number, the used method and the alarms.

### TECHNICAL SPECIFICATIONS

**Optical group:** Made up by a light source with automatic adjustment of its intensity, a diffusion lens, a diaphragm to control the light path, a reaction cell and a high sensitivity revealing element which detects all O.D.variations in the center of the cuvette.

**Thermostatic Block:** Thermostated device at  $37 \pm 0,1^\circ\text{C}$  that includes the optical system and the magnetic stirrer motor. It allows to thermostate up to 15 cuvettes and one vial of reagent.

**Numeric keyboard:** Easy to handle, to program and to perform the tests.

**Alphanumeric display:** LCD of 4 x 20 characters in which appear all the results, alarm indications, temperature control, control times and all the data of the test to be done.

**Printer:** Thermal Line Dot printer of 384 dots/line with 56 mm paper width.

**Storage capacity:** Up to 10 methods.

**Program capacity:** Incubation times, waiting time, security time, calibration curves up to nine points for the methods that require it.

**Calculations:** It automatically calculates results related to the stored calibration curves with the corresponding units, activity and International Normalized Ratio. (INR)

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