

TLD Materials Specifications

Handling and thermal treatment

Consistent, well-controlled and repeatable procedures are the key to successful TLD. Variations in annealing temperature will affect dosimeter sensitivity, for example. The following guidelines are advisable to optimize the reproducibility of bare dosimeters.

Handling

Vacuum tweezers should always be used. (Avoid mechanical tweezers or fingers). Small scratches, loss of mass or foreign deposits affect light emission).

Cleaning

Rinse the dosimeters in analytical grade anhydrous methyl alcohol between normal uses. (Do not soak). Dry by leaving to evaporate for at least one hour. Anneal once before actual use, accurately following the established procedure. The anneal will also assist in removing any residual methyl alcohol.

Annealing

For annealing temperatures up to 400 °C, the containers should be made from high temperature stainless steel or aluminum.

oxidized aluminum, preferably thin to assist rapid cooling following annealing. (Do not use non-oxidized aluminum). The use of a dedicated annealing oven reduces the risk of contamination by foreign material. Place the annealing containers on open oven racks with air space all round to avoid inconsistent heat gradients. (Do not stack containers or allow them to touch the oven walls).

CAUTION - Sensitivity to Ultraviolet Light

Calcium Fluoride Dysprosium (TLD-200), Aluminum Oxide (TLD-500) and Calcium Sulfate Dysprosium (TLD-900) are extremely sensitive to UV light.

These materials should be handled and used in the absence of UV light and stored in opaque containers. Calcium Fluoride Manganese (TLD-400) is moderately UV light sensitive.

Limiting temperatures

Significance



limit for LiF/Mg/Cu/P materials
limit for PTFE encapsulation
limit for Kapton encapsulation

Temperature

- 240 °C
- 300 °C
- 400 °C



Product Overview



© 2004 Thermo Electron Corporation. All rights reserved. Question everything. Answer. *Detailed Measure Control* are trademarks of Thermo Electron Corporation. *LiTLD* Materials (2005-4-V1.0)

This specification sheet is for informational purposes only and is subject to change without notice. Thermo makes no warranties, expressed or implied, in this product summary.

Rest of Europe:

UK:
Bath Road
Bentham, Reading RG7 5RR
D31065, Erlangen
Germany
+44 (0) 118 971 221
+44 (0) 118 971 285 fax
+44 (0) 19131 93205 fax

Other Countries Worldwide:

Vittrastrasse 5
D-42293 Wermelskirchen
Germany
+49 (0) 21 96 206 0
+49 (0) 21 96 728 24 75 fax

www.thermo.com/rmp

Thermo
ELECTRON CORPORATION

Analyze • Detect • Measure • Control™

Thermo
ELECTRON CORPORATION

System overview

Single element dosimeters and assemblies are widely used in many installations and processed using Harshaw TLD Systems. These systems include a range of Readers and Irradiators with compatible software to implement calibration, radiation evaluation and management, dose algorithms, health physics record keeping, finder database, glow curve analysis and chain of custody monitoring.

Card dosimeters

2, 3 or 4 TLD elements are assembled into rigid aluminum cards and mounted within shielded filter-holders.

EXT-RAD dosimeters

Featuring cold sterilization with efficient handling and processing the System comprises:

- Barcoded 1- or 2-element chipstrate dosimeters
- Adjustable, reusable finger rings with elements in sealed pouches
- Barcoded carrier cards for readout after exposure
- Ring sealing and chip extraction
- Carrier cards for readout after exposure

DXT-RAD dosimeters

Featuring hot or cold sterilization with fast readout, the System comprises:

- TLD-100/TLD-700 disk dosimeters with novel and permanent individual barcodes
- Sealed disposable finger rings - sterilized hot or cold
- Ring sealing and chip extraction

Material features

- Available in the form of powders and accurately machined, optically transparent disks, rods, chips and cubes
- Accurate for X, gamma, beta, electron and neutron radiations
- According to choice of material
- Simulate 'point detector' in medical physics applications
- Reusable hundreds of times
- Independent of dose rate up to 1000 MGy/s
- Long-term response retention
- Powder form only for TLD-900 (CaSO₄:Dy)

Special features of TLD-100 (LiF:MG,Tl)

- Nearly tissue-equivalent
- $\pm 15\%$ sample-to-sample uniformity
- Repeatability to within 2% or better

Filter-holder features

- Gasket-sealed to exclude dirt and moisture.
- Polarized to eliminate incorrect card insertion.
- Equipped with tamper-evident seals.
- Provided with visual indication of card barcode ID through window.
- Color-coded, per customer specification; and barcoded if so specified.

Configuration	Purpose	Filter	Radiation Fields & Mixtures	Measurement Range
1	Deep dose	AlSi: 600 mg/cm ² or 1000 mg/cm ²	LiF: Mg,Ti or LiF: Mg,Cu,P	10 μGy-20 Gy 1 μGy-20 Gy
3	Skin dose	Energy Discriminator	LiF: Mg,Ti Mylar: 333 mg/cm ² Mylar: 17 mg/cm ²	7 μGy-20 Gy — 10 μGy-20 Gy
4	Lens-of-eye dose Neutron discriminator	AlSi: 300 mg/cm ² or 600 mg/cm ² AlSi: 500 mg/cm ² or 1000 mg/cm ²	LiF: Mg,Ti or LiF: Mg,Cu,P	50 μGy-20 Gy 5 μGy-20 Gy —

Card dosimeters



DXT-RAD dosimeters



EXT-RAD dosimeters



TLD Materials

Features and Technical Specifications