

Your requirements for cuffs vary by care area and patient need. Statcorp Medical cuffs may be ordered with connectors that make them compatible with most brands of automated blood pressure monitors and wall or portable manometers.

Durable, comfortable, accurate, and economical

Three styles are available:

- UltraCheck Reusable
- Manufactured without PVC/ DEHP
 - Anti-microbial protection built into the fabric*
 - Fluid barrier

- SoftCheck
- Soft fabric prevents perspiration under cuff

- SoftCheck Vinyl
- Soft vinyl is comfortable yet can be wiped off
 - Now available in Yellow Vinyl for Isolation areas

Features


- Meets AHA/ AAMI Standards
- Color coded
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Statcorp adapters allow for compatibility with monitors worldwide





UltraCheck® & SoftCheck®

Fittings



Manual bulbs



Yellow Vinyl



For Isolation Patients

| Description | Single tube | Dual tube | Limb circumference |
|--------------------------------|-------------|-----------|--------------------|
| UltraCheck Infant | US0814-- | UD0814-- | 08 – 14 cm |
| UltraCheck Child | US1320-- | UD1320-- | 13 – 20 cm |
| UltraCheck Small Adult | US1826-- | UD1826-- | 18 – 26 cm |
| UltraCheck Adult | US2635-- | UD2635-- | 26 – 35 cm |
| UltraCheck Adult Long | US2938-- | UD2938-- | 29 – 38 cm |
| UltraCheck Large Adult | US3242-- | UD3242-- | 32 – 42 cm |
| UltraCheck Large Adult Long | US3544-- | UD3544-- | 35 – 44 cm |
| UltraCheck Thigh | US4250-- | UD4250-- | 42 – 50 cm |
| SoftCheck Infant | ST0814-- | DT0814-- | 08 – 14 cm |
| SoftCheck Child | ST1320-- | DT1320-- | 13 – 20 cm |
| SoftCheck Small Adult | ST1826-- | DT1826-- | 18 – 26 cm |
| SoftCheck Adult | ST2635-- | DT2635-- | 26 – 35 cm |
| SoftCheck Adult Long | ST2938-- | DT2938-- | 29 – 38 cm |
| SoftCheck Large Adult | ST3242-- | DT3242-- | 32 – 42 cm |
| SoftCheck Large Adult Long | ST3544-- | DT3544-- | 35 – 44 cm |
| SoftCheck Thigh | ST4250-- | DT4250-- | 42 – 50 cm |
| SoftCheck Vinyl Infant | VST0814-- | VDT0814-- | 08 – 14 cm |
| SoftCheck Vinyl Child | VST1320-- | VDT1320-- | 13 – 20 cm |
| SoftCheck Vinyl Small Adult | VST1826-- | VDT1826-- | 18 – 26 cm |
| SoftCheck Vinyl Adult | VST2635-- | VDT2635-- | 26 – 35 cm |
| SoftCheck Vinyl Adult Long | VST2938-- | VDT2938-- | 29 – 38 cm |
| SoftCheck Vinyl Large Adult | VST3242-- | VDT3242-- | 32 – 42 cm |
| SoftCheck Vinyl Lrg Adult Long | VST3544-- | VDT3544-- | 35 – 44 cm |
| SoftCheck Vinyl Thigh | VST4250-- | VDT4250-- | 42 – 50 cm |
| SoftCheck Yellow Vinyl Small | YVST1826-- | VDT1826-- | 18 – 26 cm |
| SoftCheck Yellow Vinyl Adult | YVST2635-- | VDT2635-- | 26 – 35 cm |
| SoftCheck Yellow Vinyl Large | YVST3242-- | VDT3242-- | 32 – 42 cm |

--To order cuffs for a specific monitor, add the appropriate suffix to the product number. For Critikon Dinamap®, add DM, for Datascope®, Spacelabs®, and Philips®, add HP, for GE Healthcare®, add MQ, and for Datex- Ohmeda®, add MQMF.

*The mechanism of action of this agent does not contribute to the development of drug-resistant microbial strains. Reduction in colonization or microbial growth on the device has not been shown to correlate with a reduction in infections in patients. Clinical studies to evaluate reduction in infection have not been performed.

Statcorp Medical
 an OSI Optoelectronics company
 14476 Duval Place West, Suite 303
 Jacksonville, FL 32218
 800-992-0014, 904-786-5113
 www.statcorpmedical.com



ISO 13485
 Certificate Number FM 40458

Dealer Imprint

Statcorp Medical reserves the right to make changes to this sheet and the product at any time without notice.

Velcro® is a registered trademark of VELCRO USA, Inc. All other product names mentioned herein may be trademarks and/or registered trademarks of their respective owners.